Alpinglow Summit Bulb

Specification

Power: 54w
Spectrum: 660nm850nm
Dimensions: 12x12.5cm
Material: Aluminum Alloy
LED Quantity: 18pcs 3W Led
Beam angle: 30 degree
Power Consumption: 30w+-5w
Power Input: AC 85-265V(50/60HZ)
Base type: fits both E26&E27
Lifespan: 50,000hrs

PAR38
E26 Socket

Using this product

Simply screw the device into any suitable E26/E27 socket, like a lamp or cable, and turn it on as you would any bulb.

When powered on, the 660nm LEDs emit a bright red light that is clearly visible. In contrast, the 850nm LEDs do not appear to be lit, as near-infrared light falls outside the visible spectrum and cannot be seen by the human eye.

Basic application:

- Distance: Position the light 0 to 30 cm from the target area.
- Session Duration: Expose each area for 5 to 25 minutes, depending on the distance. Repeat on other areas as needed.
- Frequency: Use 2 to 14 times per week. Listen to your body and adjust time and frequency based on how you feel.
- Note: For long-term benefits, continue using the therapy several times per week.

Thank you for investing in your health and well-being with Alpinglow Summit Bulb.

www.alpinglow.ch

[alpinglow_switzerland](https://www.instagram.com/alpinglow_switzerland)

[Alpinglow](https://www.facebook.com/Alpinglow)

info@alpinglow.ch

+41 79 674 71 07

Crans Montana, Switzerland

www.alpinglow.ch

Disclaimer:

Our guide and red light therapy panels are intended for general wellness purposes only. They are not a substitute for professional medical advice, diagnosis, or treatment. Individual results may vary. Always consult a qualified healthcare professional before beginning any new therapy or wellness practice.

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